PERSONAL HEALTH SCREENING SYSTEM KOLIBRI®

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MD KOLIBRI is Active Medical Device of **Class IIa**, according to the ANNEX VIII, Chapter I DURATION OF USE: «Transient», ACTIVE DEVICES: «Active device intended for diagnosis and monitoring» **ANNEX VIII, Chapter III 6.2 Rule 10** and SaMD KOLIBRI **6.3. Rule 11** of the REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL

MD KOLIBRI is Active Medical Device of **Class IIa**, according to the **Annex IX Sect. I clause 1.6, 3.2 Rule 10** of the Directive 2007/47/EC of the European Parliament and the Council amending Council Directive 93/42/EEC concerning medical devices SaMD KOLIBRI Class IIa according to the **ANNEX IX Sect. II clause 2.3** – of the Directive 2007/47/EC of the European Parliament and the Council amending Council Directive 93/42/EEC concerning medical devices

SaMD KOLIBRI **Class A** according to the IEC 62304:2006/AMD 1:2015 (IEC 62304:20015 ed.1.1) Medical device software - Software life cycle processes -Amendment 1

TR 753-UA Annex 2 Clause 18 - Active medical device, designated for diagnostics, are classified up to class IIa.

According to CISPR 11: 2015 + AMD1: 2016 + AMD2: 2019 CSV – Class B

GMDN 46565

UMDNS 12461

Note: SaMD - Software as a medical device Requirements for FDA, EU

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General Health Parameters (7)

| Name | Description |
|--------------|---|
| HR | Heart Rate is the speed of the heartbeat; measured by the number of contractions (beats) of the heart per minute (bpm). The heart rate can vary according to the body's physical needs, including the need to absorb oxygen and excrete carbon dioxide. |
| Health Index | Health index is a number based on scientific calculations from 0 (low index) to 10 (high index), when 0 (zero) is an excellent result. It grows or decreases in real-time and depends on changes in the parameters of your body, emotional well-being and lifestyle. |
| Stress Index | The stress index reflects the degree of stress on the body. Stress is a response of the human body to an overstrain, negative emotions, or simply to monotonous bustle. Stress in small quantities is needed by everyone, as it makes you think, to seek a way out of the problem. But severe stress affects health. Immunity decreases and a number of diseases develop (cardiovascular, gastrointestinal, etc.). It reflects the degree of centralization of heart rhythm control and characterizes the activity of the sympathetic department of the autonomic nervous system (the degree to which the activity of central regulation mechanisms prevails over autonomous ones). |
| BMI | Body Mass Index (kg/m^2) is a measure of body fat based on height and weight that applies to adult men and women. The BMI is a convenient rule of thumb used to broadly categorize a person as underweight, normal weight, overweight, or obese based on tissue mass (muscle, fat, and bone) and height. |
| BMR | Basal metabolic rate (kcal/day) is a minimum number of calories required for (basal) life-sustaining functions of a body at rest (like breathing, circulation, nutrient processing, cell production) daily, without the influence of physical activity or other impacts. The value is stable within a period of time when there are no great changes in the weight of a person |
| EE | Energy Expenditure is the number of calories that a body burns at a moment of time. This value fluctuates during a day depending on physical activity |
| MAP | The mean arterial pressure (MAP) is an average blood pressure in an individual during a single cardiac cycle. It is used by doctors to check whether there's enough blood flow, resistance, and pressure to supply blood to all major organs. Doctors usually consider anything between 60 and 100 mmHg to be normal. MAP in this range indicates that there's enough consistent pressure in the arteries to deliver blood throughout the body. |

CBC and CHEMISTRY (13)

| Name | Description |
|------|--|
| HGB | Haemoglobin (HGB) is the iron-containing oxygen-transport metalloprotein in the red blood cells (erythrocytes). Haemoglobin in blood carries oxygen from the lungs or gills to the rest of the body (i.e. the tissues). It also transports carbon dioxide out of your cells and back to your lungs to be exhaled. |
| HCT | The haematocrit (HCT) is the volume percentage (vol%) of red blood cells (RBC) in blood. Its value depends on the number and size of red blood cells. Because the purpose of red blood cells is to transfer oxygen from the lungs to body tissues, a blood sample's haematocrit—the red blood cell volume percentage—can become a point of reference of its capability of delivering oxygen. |

| | |
|------|---|
| | Haematocrit levels that are too high or too low can indicate a blood disorder, dehydration, or other medical conditions. |
| PLT | Platelets, also called thrombocytes are a component of blood whose function (along with the coagulation factors) is to react to bleeding from blood vessel injury by clumping, thereby initiating a blood clot. The number of platelets varies across individuals. The normal physiologic range is 200000 to 500000 per microliter of blood. |
| RBC | Red blood cells (RBCs), also referred to as erythrocytes, are the most common type of blood cell and the principal means of delivering oxygen to the body tissues—via blood flow through the circulatory system. RBCs take up oxygen in the lungs and release it into tissues while squeezing through the body's capillaries. |
| WBC | White blood cells (WBCs), also called leukocytes, are the cells of the immune system that are involved in protecting the body against both infectious disease and foreign invaders. All white blood cells are produced and derived from multipotent cells in the bone marrow known as hematopoietic stem cells. Five main types are neutrophils, eosinophils (acidophiles), basophils, lymphocytes, and monocytes. |
| MCHC | The mean corpuscular haemoglobin concentration (MCHC) is a measure of the concentration of haemoglobin in a given volume of packed red blood cell. |
| MCV | Mean corpuscular volume (MCV) is the average volume of red cells. MCV allows classification as either a microcytic anaemia, normocytic anaemia or macrocytic anaemia. |
| MCH | Mean corpuscular haemoglobin level refers to the average mass of haemoglobin (Hb) per red blood cell (RBC) in a sample of blood. MCH value is diminished in hypochromic anaemias. |
| BCI | Blood carbohydrate index (BCI) is an indicator that characterizes the patient's glucose tolerance for patients with type 2 diabetes. BCI level directly correlates with blood sugar. If the BCI is higher than the reference values (Normal ranges: 49.2-70.8), then the patient's blood glucose is considered above the normal value. And consequently, when BCI is lower than the reference values, then the patient's blood glucose is below the norm. |
| K | Potassium (K ⁺) is one of the electrolytes in the body's, which are minerals that carry an electric charge when dissolved in body fluids such as blood. The normal range of serum potassium is 3.5 - 5.0 mmol/L (3.5 and 5.0 mEq/L). |
| Na | Sodium (Na) is both an electrolyte and mineral. It helps keep the water (the amount of fluid inside and outside the body's cells) and electrolyte balance of the body. Sodium is also important in how nerves and muscles work. Most of the sodium in the body (about 85%) is found in blood and lymph fluid. |
| Ca | Calcium is the most common mineral in the body and one of the most important. The body needs it to build and fix bones and teeth, help nerves work, make muscles squeeze together, help blood clot, and help the heart to work. Almost all of the calcium in the body is stored in bone. |

LIPID METABOLISM (4)

| Name | Description |
|------|-------------|
|------|-------------|

| | |
|-------|---|
| TG | Triglycerides are the main constituents of body fat in humans and other vertebrates, as well as vegetable fat. They are also present in the blood to enable the bidirectional transference of adipose fat and blood glucose from the liver and are a major component of human skin oils. |
| LDL-C | Low-density lipoprotein (LDL) is one of the five major groups of lipoprotein which transport all fat molecules around the body in the extracellular water. These groups, from least dense, compared to surrounding water (largest particles) to most dense (smallest particles), are chylomicrons. LDL delivers fat molecules to the cells and can drive the progression of atherosclerosis if they become oxidized within the walls of arteries. Often, it is called "bad cholesterol", because a high LDL level leads to a build-up of cholesterol in the arteries. But it is an essential transport system for lipids the human body needs to survive. Even "small" LDL is necessary to conduct nutrients to vessels that "large" LDL can't reach. |
| HDL-C | HDL stands for high-density lipoproteins. It is sometimes called the "good" cholesterol because it carries cholesterol from other parts of a body back to the liver. |
| CHOL | Cholesterol is an organic molecule, a sterol (or modified steroid), a type of lipid. Cholesterol is biosynthesized by all animal cells and is an essential structural component of animal cell membranes. Cholesterol also serves as a precursor for the biosynthesis of steroid hormones, bile acid and vitamin D. Hepatic cells typically produce the greatest amounts of cholesterol. |

HEART ACTIVITY (11)

| Name | Description |
|-------------|---|
| RR interval | The time between beats is measured in milliseconds (ms) and is called an "R-R interval" or "inter-beat interval (IBI)." While heart rate focuses on the average beats per minute, heart rate variability (HRV) measures the specific changes in time (or variability) between successive heartbeats. |
| P-wave | The P wave represents the depolarization of the left and right atrium and also corresponds to atrial contraction and is usually 0.08 to 0.10 seconds (80-100 ms) in duration. |
| PR segment | The PR segment is the flat line between the end of the P-wave and the start of the QRS complex. The PR segment reflects the time delay between atrial and ventricular activation. The PR segment also serves as the baseline (reference line or isoelectric line) of the ECG curve. |
| PR interval | The PR interval is the time from the onset of the P wave to the start of the QRS complex. It reflects conduction through the AV node. The normal PR interval is between 120 – 200 ms (0.12-0.20s). |
| QRS complex | The QRS complex is a name for the combination of three of the graphical deflections seen on a typical electrocardiogram. It is usually the central and most visually obvious part of the tracing; in other words, it's the main spike seen on an ECG line. It corresponds to the depolarization of the right and left ventricles of the human heart and contraction of the large ventricular muscles. In adults, the QRS complex normally lasts 0.06–0.10 s; in children and during physical activity, it may be shorter. |
| QT interval | The QT interval is calculated as the time from the start of the Q wave to the end of the T wave and approximates to the time taken from when the cardiac ventricles start to contract |

| | |
|----------------------|---|
| | to when they finish relaxing. An abnormally long or abnormally short QT interval is associated with an increased risk of developing abnormal heart rhythms and sudden cardiac death. |
| ST segment | The ST segment represents the isoelectric period when the ventricles are in between depolarization and repolarization. The ST segment connects the QRS complex and the T wave and has a duration of 0.005 to 0.150 sec (5 to 150 ms). The typical ST segment duration is usually around 0.08 sec (80 ms). It should be essentially level with the PR and TP segments. The most important cause of ST segment abnormality (elevation or depression) is myocardial ischaemia or infarction. |
| ST interval | In an electrocardiogram, the interval that represents ventricular repolarization. |
| T-wave | In electrocardiography, the T wave represents the repolarization of the ventricles. The T wave contains more information than the QT interval. The T wave can be described by its symmetry, skewness, slope of ascending and descending limbs, amplitude and subintervals like the Tpeak–Tend interval. |
| Morphology of T-wave | T-wave morphology analysis (TMA) quantifies irregularities of ventricular repolarization. TMA is useful for risk stratification of patients with myocardial infarction (MI). |
| QTc interval | The corrected QT interval (QTc) estimates the QT interval at a standard heart rate of 60 bpm. This allows comparison of QT values over time at different heart rates and improves detection of patients at increased risk of arrhythmias. |

| Heart Rate Variability: (30) | |
|------------------------------|--|
| Name | Description |
| ARP index | Index of adequacy of the regulatory processes |
| AR index | Autonomic rate index |
| AB index | Autonomic balance index |
| SER index | Normalized index of summary effect of regulation |
| ASER index | Conditional index of activity of the sympathetic element of regulation |
| FCS index | The most probable level of functioning of cardiovascular system |
| MARI | Maximum amplitude of regulatory influences |
| TINN | Triangular interpolation of the normal-to-normal interval of a histogram |
| HRVti | HRV Triangular Index. It demonstrates the general variability of the heart rhythm and is directly proportional to parasympathetic activity. |
| sdHR | Deviation of heart rate |
| RMSSD | The square root of the mean of the sum of the squares of differences between adjacent NN intervals. It allows estimating the short-term components of HRV. RMSSD characterizes the activity of parasympathetic autonomic regulation. |
| pNN50 | The proportion derived by dividing NN50 by the total number of NN intervals, when NN50 is the number of interval differences of successive NN cardiac intervals greater than 50 ms, %. The degree of prevalence of the parasympathetic regulation over the sympathetic (relative value). |
| SDNN | Standard deviation of the "normal-to-normal" intervals. SDNN characterizes the total effect of the sympathetic and parasympathetic systems on the regulation of cardiac activity (reduced |

| | |
|---------|--|
| | – predominance of sympathetic regulation; elevated – predominance of autonomic regulation). |
| MeanRR | Mean RR-interval |
| MinRR | Minimum RR-interval |
| MaxRR | Maximum RR-interval |
| SampEn | Sample entropy measures the regularity and complexity of a time series |
| Ellipse | Area of the ellipse which represents total HRV (Poincaré plot) |
| SD1 | Poincaré plot standard deviation perpendicular the line of identity |
| SD2 | Poincaré plot standard deviation along the line of identity |
| DFA1 | Detrended fluctuation analysis, which describes short-term fluctuations (α_1) |
| DFA2 | Detrended fluctuation analysis, which describes long-term fluctuations (α_2) |
| nLF | Normalized power of the low-frequency band (0.04–0.15 Hz), also called Relative power of the low-frequency band (%) |
| nHF | Normalized power of the high-frequency band (0.15–0.4 Hz), also called Relative power of the high-frequency band (%) |
| LF/HF | Ratio of LF-to-HF power (%) |
| VLF | The absolute power of the very-low-frequency (VLF) band (0.0033–0.04 Hz) reflects even slower modulations of heart rate (ms^2). The VLF band may represent the influence of the peripheral vasomotor and renin-angiotensin systems |
| LF | The absolute power of the low-frequency (LF) band (0.04–0.15 Hz) |
| HF | The absolute power of the high-frequency (HF) band (0.15–0.4 Hz) |
| TP | Total power (TP) is the sum of the energy in the VLF, LF, and HF bands for short-term recordings (less than 24h) |
| RSA | Respiratory sinus arrhythmia (RSA) is heart rate variability in synchrony with respiration, by which the R-R interval on an ECG is shortened during inspiration and prolonged during expiration. Although RSA has been used as an index of cardiac vagal function. |

| RISK OF PATHOLOGIES (14) based on ECG findings + HRV diagnostics | | displayed only if identified |
|---|--|------------------------------|
| ICD-10 | Description | |
| I21 | ST elevation (STEMI) and non-ST elevation (NSTEMI) myocardial infarction | |
| R94.31 | Abnormal electrocardiogram | |
| I44.0 | Atrioventricular block, first degree | |
| I44.5 | Left posterior fascicular block | |
| I44.7 | Left bundle-branch block, unspecified | |
| I45.6 | Pre-excitation syndrome | |
| I45.9 | Conduction disorder, unspecified | |
| I44.4 | Left anterior fascicular block | |
| I45.1 | Other and unspecified right bundle-branch block | |
| R00 | Tachycardia, unspecified | |
| R00.1 | Bradycardia, unspecified | |
| I24 | Other acute ischemic heart disease | |
| I45.81 | Long QT syndrome | |
| I48 | Atrial fibrillation and flutter | |

| ECG FINDINGS (43) based on HA | | displayed only if identified | |
|----------------------------------|------|--|---|
| Method and corresp. Code | | ICD10 | Description |
| Simplified Minnesota Code | 1.11 | I21 | ST elevation (STEMI) and non-ST elevation (NSTEMI) myocardial infarction |
| Simplified Novacode Code | 5.1 | | Q wave MI; major Q waves with or without ST-T abnormalities |
| Simplified Minnesota Code | 1.12 | | ST elevation (STEMI) and non-ST elevation (NSTEMI) myocardial infarction |
| Simplified Novacode Code | 5.2 | | Q wave MI; moderate Q waves with ST-T abnormalities |
| Simplified Minnesota Code | 1.21 | | ST elevation (STEMI) and non-ST elevation (NSTEMI) myocardial infarction |
| Simplified Novacode Code | 5.3 | | Possible Q wave MI; moderate Q waves without ST-T abnormalities |
| Simplified Minnesota Code | 1.22 | | ST elevation (STEMI) and non-ST elevation (NSTEMI) myocardial infarction |
| Simplified Novacode Code | 5.4 | | Possible Q wave MI; minor Q waves with ST-T abnormalities |
| Simplified Minnesota Code | 4.3 | R94.31 | Minor ST-T abnormalities |
| Simplified Minnesota Code | 4.4 | | Minor ST-T abnormalities |
| Simplified Novacode Code | 5.5 | | ST abnormalities without Q waves |
| Simplified Minnesota Code | 5.1 | | T wave abnormalities without Q waves / Possible Q wave MI; minor Q waves with ST-T abnormalities / Q wave MI; moderate Q waves with ST-T abnormalities |
| Simplified Novacode Code | 5.6 | | T wave abnormalities without Q waves |
| Simplified Minnesota Code | 5.2 | | Left ventricular hypertrophy with ST-T / T wave abnormalities without Q waves / Possible Q wave MI; minor Q waves with ST-T abnormalities / Q wave MI; moderate Q waves with ST-T abnormalities |
| Simplified Novacode Code | 5.7 | | Minor Q waves without ST - T abnormalities |
| Simplified Minnesota Code | 5.3 | | Minor ST-T abnormalities |
| Simplified Minnesota Code | 5.4 | | Minor ST-T abnormalities |
| Simplified Novacode Code | 5.8 | | Minor ST-T abnormalities |
| Simplified Novacode Code | 1.41 | I47.1 | Supraventricular tachycardia, rate < 130 cpm |
| Simplified Novacode Code | 1.42 | | Supraventricular tachycardia, rate ≥ 130 cpm |
| Simplified Minnesota Code | 6.3 | I44.0 | First-degree AV block (AVB1) |
| Simplified Novacode Code | 2.1 | | Atrioventricular block, first degree |
| Simplified Minnesota Code | 7.7 | I44.4 | Left anterior fascicular block (LAFB) |
| | 7.11 | | Left bundle branch block without ECG evidence of myocardial infarction (MI)/ QRS duration ≥ 120 ms |
| Simplified Novacode Code | 3.1 | | Left bundle branch block without ECG evidence of myocardial infarction |
| Simplified Novacode Code | 3.11 | | Left bundle branch block with possible myocardial infarction (MI) |
| Simplified Minnesota Code | 7.21 | I45.1 | Right bundle branch block without ECG evidence of MI/ Right bundle branch block with possible MI/QRS duration ≥ 120 ms |
| Simplified Novacode Code | 3.2 | | Right bundle branch block without ECG evidence of MI |
| Simplified Minnesota Code | 6.6 | I45.9 | Intermittent aberrant atrioventricular conduction |
| Simplified Minnesota Code | 7.4 | | Indeterminate ventricular conduction delay without ECG evidence of MI/ QRS duration ≥ 120 ms |
| Simplified Novacode Code | 3.1 | | Indeterminate ventricular conduction delay without ECG evidence of MI |
| Simplified Novacode Code | 3.31 | | Indeterminate ventricular conduction delay with possible MI |
| Simplified Novacode Code | 3.41 | | Borderline delay of right ventricular excitation |
| Simplified Novacode Code | 3.42 | | Borderline delay of left ventricular excitation |
| Simplified Novacode Code | 4.11 | | I45.81 |
| Simplified Novacode Code | 4.12 | Significant prolongation of ventricular repolarization | |
| Simplified Minnesota Code | 3.1 | I51.7 | Left ventricular hypertrophy without ST-T |
| Simplified Minnesota Code | 4.11 | | Left ventricular hypertrophy with ST-T/ST abnormalities without Q waves |
| Simplified Minnesota Code | 4.12 | | Left ventricular hypertrophy with ST-T/ST abnormalities without Q waves |

| | | | |
|---------------------------|------|--------|--|
| Simplified Minnesota Code | 4.2 | | Left ventricular hypertrophy with ST-T/ST abnormalities without Q waves |
| Simplified Minnesota Code | 6.41 | I44.7 | QRS duration \geq 120 ms/Ventricular preexcitation pattern (WPW) |
| Simplified Minnesota Code | 6.5 | I45.6 | Short P-R interval |
| Simplified Minnesota Code | 7.6 | I44.5 | Left posterior fascicular block (LPFB) / Borderline delay of left ventricular excitation |
| Simplified Minnesota Code | 8.7 | R00 | Sinus Tachycardia (ST) |
| Simplified Minnesota Code | 8.8 | R00.1 | Sinus Bradycardia (SB) |
| Simplified Minnesota Code | 9.2 | I24 | Infarction / Ischemia |
| Simplified Minnesota Code | 9.1 | R94.31 | Low QRS amplitude |
| Simplified Minnesota Code | 9.5 | | T-wave amplitude > 12 mm |

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| | |
|--|--|
| Mains voltage (via USB port) | U = 5±1 V DC I ≤ 800 mA |
| Battery (Lithium Polymer Cell) | LP502060 3.7V 510mAh, installed on the PCB (see the details below) |
| Consumption | max. 2,5 VA |
| MD Class | II A |
| Basic safety/EMC | IEC 60601-1:2005+AMD1:2012 (ed.3)/IEC 60601-1-2 4th Ed |
| Class of protection against penetration of water or solid particles (IP) | IP 20 |
| Readiness time counted from switching on | max. 10 sec |
| Net Weight | 55 ±5 g |
| Dimensions | 82x56x10 mm |
| Length of cables with electrodes | 0.70-0.85 m (set of cables is optional) |
| Time of examination | 5 min |
| Applied methods | Definition of human health parameters on the base of HA measured at designated points on a body (on the skin) and anthropometric data of the patient. |
| Software | KOLIBRI APP (installed on a mobile gadget), KOLIBRI SW (cloud based), the hardware has internal software, preinstalled on a chip |
| Ambient temperature for storage and transportation | From -5°C (without relative humidity control) till + 45°C (with relative humidity control) Class 7K2, as described in the IEC/TR 60721-4-7:2001 |
| Ambient temperature and humidity during operation of MD | The device is to be used indoor. Environmental conditions should not induce transient conditions in a person such as excessive sweating or shivering from cold. Therefore, a |

| | |
|--|---|
| | comfortable temperature and humidity in the room must be ensured. |
| Atmospheric pressure (operation) | 650-765 mm Hg (87 -102 kPa). |
| Humidity | 15%-93%, without condensation |
| Temperature range | +20 °C to +27 °C |
| Humidity of the surface of skin during testing with MS KOLIBRI | 90-95% - the skin in places of contact must be moisturized before taking the measurement (ECG gel is preferred; water is possible) |
| Precision pressure sensor with altimetry | Pressure absolute accuracy (Test conditions: 50 to 110 kPa over -10 °C to 70 °C) - ± 0.4 kPa |
| | Pressure relative accuracy (Relative accuracy during changing temperature between -10 °C to 50 °C at any constant pressure between 50 kPa to 110 kPa) |
| The temperature of all surfaces the medical device: | The temperature of all surfaces of the medical device does not exceed 40 ± 0.5 °C |

Minimum requirements to a smartphone or tablet, which are enough to use with MS KOLIBRI

| Minimum hardware requirements | | Minimum software requirements | |
|-------------------------------|---|---|-------------------------|
| Display | min. 5"-55" definition HD+, FHD+, QHD+, WQHD+ or more | Operating system | Android OS 6 or higher. |
| Mobile gadget support | LTE, 3G or high, Wi-Fi, BLUETHOOH 3.0 or high | Web browser | Any web browser |
| Antivirus software | | Any antivirus software must be set up in such a way when the software KOLIBRI is not blocked. | |

As the KOLIBRI has a battery for wireless mode of operation, the following technical features shall be considered:

| Battery (Lithium Polymer Cell) LP502060 3.7V 510-530 mAh, installed on the PCB | |
|---|---|
| Charge limited voltage/ Nominal voltage | 4.2V/3.7V (<i>the cell is protected by an electronic circuit that won't allow it to overcharge nor over-discharge under use</i>) |
| Rated capacity | 500mAh min, 530mAh typ. |
| Wat-Hour Rating | 1.87 Wh |
| Max. Operating Voltage Range | 2.7V to 4.20V |
| Max. Charge Voltage | 4.2V ± 50 mV |
| Max. Charge Current | 500mA |
| Discharge Current | 500mA |
| Discharge Cut Off | 2.75V |
| Cycle life (0.5C/0.5C) 23 ± 5 °C | 500 times (<i>approximately 1 year of operation</i>) <i>One cycle refers to one charge period and then one discharge period. So, cycle life means that after 500 times of discharge-charge, the cell capacity will be reduced to 80% of the rated one.</i> |
| Cell Protection | |
| Overcharge Detection | 4.275 ± 50 mV (0.7 to 1.3sec. delay, release 4.275V ± 50 mV) |
| Over-discharge Detection | 2.75V ± 50 mV (14 to 26msec. delay, resume 2.50V ± 50 mV) |
| Overcurrent Detection | 2A to 4.5A (8 to 16msec. delay) |
| Ambient Conditions | |
| Charge Temp. Range | 0 to +45 °C |

| | |
|--------------------------|---|
| Discharge Temp. Range | -20 to +60 °C |
| Storage Temp. Range | 1 year at -20 to +30°C >70% |
| Overcurrent Detection | 3 months at -20 to +45°C >70% 1 month at -20 to +60°C >70% |
| Humidity | 65 ±20%RH |
| Environmental and Safety | Please follow LiPo Handling and Safety Precautions for Lithium Polymer Battery. This battery meets the requirements of Battery Directives, and the battery parts are IEC 62133 & RoHS-Compliant. For more safety precautions and performance standards, please go to www.lipolbattery.com/support.html |

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
Sign in to the MS KOLIBRI is done every time after log out of the account. If you close the application and do not logout, you will be kept signed in. Therefore, the next time you run the application, sign in will not be required.

During the first login to the MS KOLIBRI, it is necessary to register your device in the system. Please press **“Sign up”** for that.

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If you understand everything and agree to the terms, please tick each of the check-boxes and click on: «ACCEPT TERMS».

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15:04 90%

format and has the right to transmit that data to another controller without hindrance from the controller to which the personal data has been provided...

For the KOLIBRI software:
The KOLIBRI software allows you to export data from the results of a patient examination onto paper and PDF files with the subsequent transfer of files via email. This enables you to realize the right to data portability.
The KOLIBRI software allows you to remotely connect the transfer of the results of your research to the doctor's cloud account and the function of displaying the test results in the doctor's account in real-time is implemented.
This allows you to realize the right to data portability.

Signature of patient

GDPR: I have been briefed on my rights and Privacy Policy based on the REGULATION (EU) 2016/679 (GDPR) regarding the processing and movement of the personal data, which are used for the purposes of non-invasive hemogram screening, provided by the software KOLIBRI. I am aware of and agree with the condition that the data containing in the database of the software KOLIBRI will not be disclosed and/or forwarded to unauthorised third parties except the following: a) to a personal data processor (a professional, who is responsible for carrying out an examination on site); b) to the respective/authorised staff of the manufacturer, whose duties includes elimination of bugs in the software KOLIBRI, software update or other assistance related to the software KOLIBRI on the request of the personal data processor; c) to governments, control agencies, regulatory authorities and others as permitted or required by law. I am aware of an ability to withdraw my consent at any time submitting a written statement in the place, where the service was provided. Therefore, I realise that my personal data will be deleted irrevocably and without a possibility of further applying, referring and analysing of such a data.

Signature of patient

INTERNATIONAL TECHNOLOGY CENTER
P.O.BOX 330102, Ras Al Khaimah, United Arab Emirates
Tel. +971-07-2041957
mob. +971-55-7702115
e-mail: itc.fzc@gmail.com
http://intertc.com/about_us

ACCEPT TERMS

14:58 61%

Ver 1908081742

Sign up KOLIBRI

Full name

Email

Password

Show password

SIGN UP

Already have an account? Sign in

Open your email box. Check the Inbox and SPAM folders. Open the letter from MS KOLIBRI and click on "**Confirm email**".

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After successful activation, you will be automatically forwarded to the menu "Medical Card"

APPLICATION KOLIBRI

17:43 100%

← Medical card

Test patient 12345

Date of birth EDIT

Gender male female

Height cm ft

Weight kg lb

Race Light-skin

Blood type Unknown

Rh factor Unknown

Diabetes Patient doesn't h..

Phone number Enter phone number

Smoking ☐

Using drugs ✓

WEB-SITE <https://kolibri.one>

16:40 39%

full name, date of birth, gender, weight, height, race and diabetes. After saving the data, fields date of birth and gender will be locked in order to avoid using of the device by unauthorized persons. To change these data, please contact the customer support: support@kolibri.one

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16:40 39%

☰

contact the customer support: support@kolibri.one

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
Select menu “FINANCE” and appropriate submenu

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The screenshot shows the 'Balance replenishment' interface. At the top, there are three columns of information: 'Available balance' (436.2 USD), 'Min. replenishment' (2 USD), and 'Max. replenishment' (2000 USD). Below this is a form with a text input field for 'Amount of replenishment (USD)' and a dropdown menu for 'Payment method' (currently set to 'PayPal'). A green 'Replenish' button is at the bottom. Red markers 1-5 point to the available balance, min/max replenishment, the input field, and the payment method dropdown respectively.

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This alternative screenshot shows the 'Amount of replenishment (USD)' field with the value '100'. Below it, the 'Payment method' section features two buttons: 'PayPal' and 'Paymentwall'. A green 'Replenish' button is at the bottom.

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KOLIBRI **UNPAID**

Invoice #152 1

| Payer name | Personal ID | Invoice date | Payment date | Payment method |
|---------------|-------------|--------------|--------------|----------------|
| Patric Burton | 1 | 01/18/2018 | — | PayPal |

Description

| | |
|-----------------------|--------------|
| Balance replenishment | 2 USD |
| Transaction amount | 0 USD |
| TO PAY | 2 USD |

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
| Date | Description | Sum |
|------------|-------------------------------------|---------|
| 01/18/2018 | Balance replenishment. Invoice #155 | + 2 USD |
| 01/18/2018 | Balance replenishment. Invoice #154 | + 2 USD |
| 01/18/2018 | Balance replenishment. Invoice #153 | + 2 USD |

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| Date | Description | Sum | Actions |
|------------|-----------------------|-----------|---------|
| 03/11/2020 | Self patient creation | - 0.5 USD | Actions |

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| No. | Subject | Invoice date | Payment date | Sum | Status | Actions |
|-----|---|--------------|--------------|---------|--------|---------|
| 11 | Balance replenishment | 10/22/2018 | — | 2 USD | UNPAID | Actions |
| 10 | Purchase of plan: Extended (365 days) | 10/22/2018 | 10/22/2018 | 357 USD | PAID | Actions |
| 9 | Purchase of plan: Starter (365 days) | 10/22/2018 | 10/22/2018 | 89 USD | PAID | Actions |
| 8 | Purchase of plan: Extended (365 days) | 10/22/2018 | 10/22/2018 | 357 USD | PAID | Actions |
| 7 | Purchasing of an activation key. Plan: Starter (365 days) | 10/22/2018 | 10/22/2018 | 89 USD | PAID | Actions |

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The screenshot displays three corporate pricing plans side-by-side. Each plan is presented in a card format with a blue header, a description, a monthly fee, a savings amount, and a list of included services. The plans are: Doctor S (a few patients), Doctor M (medium and large clinics), and Doctor L (wide private practice). The monthly fees are \$60, \$119, and \$238 respectively. Savings are 53.76 USD, 107.52 USD, and 215.04 USD. Services include tests per day, patient and staff creation costs, and free patients/staff. Each plan also shows coverage for 41 out of 41 parameters and 17 out of 17 pathologies, with a 12-month term. A 'Choose plan' button is at the bottom of each card. Red markers (1-14) and blue question marks are overlaid on the image to highlight specific UI features.

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| | | | |
|-----|-----------------|---|--|
| 10. | Низький баланс. |  <p>The screenshot shows the KOLIBRI app interface. At the top, the status bar displays the time 21:39 and battery level 49%. The app header shows 'KOLIBRI' with a back arrow and a menu icon. Below the header, a summary card displays: 'Tests 0 2 Bonuses 0', 'Failed tests today - 1', 'Ballance 0.0 USD', 'Price per test 0.2 USD', 'Subscription plan: Personal M', and 'Subscription expires in 364 days'. A white error dialog box is centered on the screen with the text 'You don't have enough money on balance' and three buttons: 'Check', 'Replenish', and 'OK'. Below the dialog is a 'Results' section with a flask icon. The bottom navigation bar shows three icons: a list, a home button, and a back arrow.</p> | <p>Поповніть баланс на своєму рахунку. Поповнивши баланс, повторіть тест.</p>  <p>The inset screenshot shows a close-up of the error message card. Red boxes highlight the 'Bonuses 0' and 'Ballance 0.0 USD' fields. The text of the error message is: 'You don't have enough money on balance'.</p> |
|-----|-----------------|---|--|

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The MD KOLIBRI

| |
|--|
| PERSONAL HEALTH SCREENING SYSTEM KOLIBRI® |
|--|

Name of MD KOLIBRI

Model

 KOLIBRI

Serial Number

is packed at the manufacturing facilities of the company SCIENTIFIC COMPANY KOLIBRI LLC in accordance with the Technical Specification

Complete set includes

| |
|--|
| |
| |
| |
| |
| |

Date of packing

| | | |
|-----------|-----------|-------------|
| <i>DD</i> | <i>MM</i> | <i>YEAR</i> |
|-----------|-----------|-------------|


Packed by

| | |
|------------------|-------------------------|
| <i>Signature</i> | <i>Name and Surname</i> |
|------------------|-------------------------|

**Manufacturer:****SCIENTIFIC COMPANY KOLIBRI LLC**

Kinnyi provulok 8a, Kharkiv, 61001, Ukraine.

Tel.: +380913011110

email: info@kolibri.oneURL: <https://kolibri.one>**Authorized representative in European Union:****Authorized representative in European Union:****HUNGARY****Medical Devices Magyarország KFT****1149 Budapest, Vezér utca 79-2****Tel: +36-20-971-93-23**email: info@mdevice.orgURL: <https://www.mdevice.org/>

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| | | |
|--|--|-------------------------------------|
| Warranty card | PERSONAL HEALTH SCREENING SYSTEM KOLIBRI® | |
| Card No. (Serial number of the MD KOLIBRI) | | |
| Warranty period* | 36 months for MD KOLIBRI | |
| Purchase Date | Seller | |
| <i>YEAR/MM/DD</i> | <i>Trade organization Name</i> | |
| | Seller address | Telephone |
| | <i>Trade organization Address</i> | <i>Trade organization telephone</i> |
| Customer | <i>[Customer Signature]</i> | <i>[Signature]</i> |
| <i>Hereby, I confirm that I have read the warranty terms and agree with them</i> | | <i>[Company stamp]</i> |

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| TEAR-OFF COUPON NO.1 | | [SERIAL NUMBER OF DEVICE] | No.1 | Filled up by service centre |
|----------------------|----------|---------------------------|-------|-----------------------------|
| PRODUCT | KOLIBRI® | | _____ | _____ |
| MODEL NAME | | | _____ | Reception date |
| PURCHASING DATE | | | _____ | Release date |
| TRADE ORGANIZATION | | | _____ | Special Marks |

| | | |
|-----------------------------|--------------------------------|---|
| | TRADE ORGANIZATION STAMP | SERVICE CENTER STAMP |
| TEAR-OFF COUPON NO.2 | [SERIAL NUMBER OF DEVICE] | No.2 Filled up by service centre |
| PRODUCT | KOLIBRI® | |
| MODEL NAME | | Reception date |
| PURCHASING DATE | | Release date |
| TRADE ORGANIZATION | | Special Marks |
| | TRADE ORGANIZATION STAMP | SERVICE CENTER STAMP |

WARRANTY CARD MAY BE PRINTED-OUT SEPERATLY FROM THE USER'S MANUAL AND INCLUDED INTO COMPLETE SET OF THE KOLIBRI TOGETHERWITH/AS A PART OF TECHNICAL PASSPORT.

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

Labels on the outer package of the MD KOLIBRI®


Labels and marking used for the MD KOLIBRI.


A number of internationally recognized symbols are placed on a label of the MD KOLIBRI®. The symbols relate to the safety requirements and standards and are briefly reviewed below:

| | | | | |
|--|--|--|----------------------|--|
| 2 | Trade name of the medical device: | KOLIBRI® | | |
| CE symbol (where xxxx is the number of Notified Body) | | | TYPE BF APPLIED PART | |
| Issue Date | [Date of manufacturing] | Serial number of MD | [Serial number] | |
| Manufactured by | SCIENTIFIC COMPANY KOLIBRI LLC Ukraine 61001 Kharkiv, Kinnyi provulok 8a, Tel.: +380913011110 Ukraine. email: info@kolibri.one URL: https://kolibri.one | Symbol for connection of USB port | | |
| Represented by (European representative) | in Europe Authorized representative in European Union: HUNGARY Medical Devices Magyarország KFT 1149 Budapest, Vezér utca 79-2 Tel: +36-20-971-93-23 email: info@mdevice.org URL: https://www.mdevice.org/ | CLASS II equipment IEC 60417-5172 | | |
| MD Class IIa | MD Class IIa | Attention label | | |
| Sign of waste electronic equipment | | Carefully read the IFU before the first use of the device | | |

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KOLIBRI®

Personal health screening system
Designed by Health Entire LTD. Assembled in EU

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KOLIBRI®

Personal health screening system
Designed by Health Entire LTD. Assembled in UA

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MD Class: IIa

According to the **Directive 2007/47/EC** of the European Parliament and of the Council of 5 September 2007 amending **Council Directive 93/42/EEC** concerning medical devices (hereinafter Directive)

KOLIBRI™ is a transient noninvasive active medical device of **Class IIa**, according to the Annex IX Sect. I clause 1.1 and clause 1.6, Sect. II clause 2.3, Sect. III clause 3.2 Rule 10 of the Directive

According to the **REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL** of 5 April 2017 on medical devices (hereinafter MDR).

KOLIBRI™ is a transient noninvasive active medical device for diagnostics of **Class IIa**, according to the Annex VIII, Sect. I clause 1.1 and clause 2.5, Sect. II clause 3.3, Sect. III clause 6.2 Rule 10 and clause 6.3 Rule 11 of the MDR.

26 COMPLIANCE WITH GDPR

Considering that on the 25th of May 2018, the EU regulation on protection of personal data entered into force, namely **REGULATION (EU) 2016/679 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 April 2016** on the protection of natural persons with regard to the processing of personal data and on the free movement of such data, and repealing Directive 95/46/EC (General Data Protection Regulation) - hereinafter **GDPR Regulation**, we, as the manufacturer of the medical equipment for screening diagnostics (**KOLIBRI™**) recommend our customers to take appropriate measures for the protection of sensitive information on the health state of their patients.

For our side, **we declare** that no information and data is transmitted from the device to third parties, as well as to us, as the manufacturer. We do not receive, store or use the information about your patients. The database is kept locally on your computer and it is your responsibility to keep it safe. Being an authorized user of the device (a holder of the software protection key), only you are able to share the information from the database of the device directly to appointed addressee. We have implemented several protective measures, which will be helpful for you in that regards, they are the following:

- ✓ Access to the database and encryption/decryption of data is carried out with a password and use SSL protocol with authorization into an email. The software and the database don't run without the password, so only the authorized user, who keeps the password to the definite account (hence to the corresponding DB) has to access to the records in the database.
- ✓ Access to the database is protected with a password (hidden) additionally and is possible only via the software **KOLIBRI™**.

We ensure, that the volume of information about a person, which is added to the database of the software is adequate, relevant and limited to what is necessary for the purposes of diagnostics and the collected data will not be further processed in a manner which is incompatible with those purposes (complies with Article 5 of the Regulation (EU) 2016/679).

We strongly advise against sending the results of testing by email to avoid any information leak or unauthorized access to the sensitive data of your patients. As you have an ability to print out the results, so you can give the survey with results of testing to your patients in a paper form, after signing the informed consent.

We recommend getting a signed consent form for processing and storage of personal data of a patient, including the data concerning his/her health, **from each of your patient BEFORE testing on the device**.

A sample of the informed consent is below. It was also included in the survey/report form and is usually printed at the beginning of the document, before the results.

INFORMED CONSENT (GDPR): *I have been briefed on my rights and Privacy Policy based on the REGULATION (EU) 2016/679 (GDPR) regarding the processing, storing and free movements of my personal data, which are used by the medical device and software KOLIBRI™ for the purposes of a non-invasive screening test.*

I am aware of and agree with the condition that the data containing in the database of the software IDEA-4RS will not be disclosed and/or forwarded to unauthorized third parties except the following:

- a) to a personal data processor (a medical professional, who is responsible for carrying out an examination and interpret the results on site);*
- b) to the respective/authorised staff of the manufacturer, whose duties includes the elimination of bugs in the software, software update or other assistance related to the software on the request of the personal data processor;*
- c) to governments, control agencies, regulatory authorities and others as permitted or required by law.*

I am aware of an ability to withdraw my consent at any time submitting a written statement in the place, where the service was provided. Thereof, I realize that my personal data will be deleted irrevocably and without a possibility of further applying, referring and analysing of such data.

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