



INSTRUCTION FOR USE (IFU)

Low-Intensity VHF-UHF Therapy Apparatus

BIOL



UA.TR.099

Version 1.0



Scientific company


KOLIBRI LLC

Provulok Kinnyi 8a, Kharkiv,
61001, Ukraine

Tel.: +380913011110

Email: info@kolibri.one

URL: <https://kolibri.cloud>



The official representative in the EU:

“ONKOCET Ltd.”

4 Kutuzovova str., 90201 Pezinok, Slovakia,

tel.: +421 (2) 44 64 09 77

Email: onkocet@onkocet.eu

URL: www.onkocet.eu

The Instructions for use, including all its sections, is protected by the copyright. Any use of materials from this manual that violates the requirements of the current copyright law without the written permission of the manufacturer is prohibited.

This provision concerns, copying, translating, micro photocopying, and placing and processing in electronic systems.

CONETNT	
Symbols and signs	4
1 General review	5
1.1 Purpose of the device	5
1.2 Indications for use	5
1.3 Contraindications	6
1.4 Side effects	7
2 Warnings and safety provisions	7
2.1 Safety signs on the device	7
2.2 Warnings and safety instructions	7
2.3 Measures to prevent damage of equipment and the device	9
3 Brief description of the device	9
3.1 Functional diagram of the device	9
4 Before using the device BIOL	11
4.1 The list of the preparatory actions	11
4.2 The operator's qualification	12
5 Operation of the device BIOL	12
5.1 A patient's location	12
5.2 An operator's location	13
5.3 Recommended operating modes	13
5.4 Switching the device on	13
5.5 Description of control means	14
5.6 Setting up the parameters of a session	15
5.7 Shutdown and storage	16
5.8 Messages on the display	16
6 Maintenance of the device BIOL	17
6.1 Cleaning of the device	17
6.2 Fuse replacement	17
6.3 Maintenance and safety check	17
6.4 Disposal and the environmental protection	18
6.5 Repairs	18
6.6 Lifetime	18
7 Completion and recommended service materials	18
7.1 Complete set	18
7.2 Recommended service materials	18
8 The device appearance	19
9 Technical descriptions	19
9.1 Classification of the device BIOL	19
9.2 Technical data	19
9.3 Storage and transportation	20
10 Warranty terms	20
11 Possible errors and malfunctions. Troubleshooting.	21
12 Checklist	21
13 Product acceptance protocol	22
14 Packing list	22
15 The warranty certificate	23
16 Libelling	24
17 The complete set of STARDEVICE®	24
18 Document history and version control	24

SYMBOLS AND SIGNS



Warning



Prohibition



Obligatory action



Sample sign for obligatory actions



Be sure to read the instruction

IP20

Degree of protection of the device against penetration of solid objects



Protective ground



Working part of type BF

ON

To turn the power on

100-240VAC

Supply voltage of the device

OFF

To turn the power off

T4AL250V

A fuse



Manufacturer

OK

Enter a value



Date of manufacture

SN

Serial number of the device



It is disposed as used electronic equipment



CLASS II
IEC 60417-5172

1 GENERAL REVIEW

1.1 INTENDED USE

Medical device BIOL can be used in clinical practice for therapy with low-intensity electromagnetic waves, clearly defined shape and frequency (100-1500MHz). Treatment is based on the restoration of damaged areas of the cell membrane and inhibition of the frequencies of pathological wave processes in the body. The medical device is designed to prevent the formation of fibrous tissues and reduce the existing fibrous formations, normalize the functioning of the immune system, treat and prevent viral diseases as part of complex therapy, reduce the duration of the postoperative and rehabilitation period, relieve pain, treat prostate pathologies.

1.2 RECOMMENDATIONS FOR USE

The device is designed to modulate the patient's immune system, treat viral diseases, reduce the duration of the postoperative rehabilitation period and relieve pain.

Individual use of the device is recommended by a doctor, who determines the duration and number of sessions, as well as monitors the course of treatment according to the patient's clinical tests.

The device is recommended for use in medical, treatment-and-prophylactic, sanatorium and outpatient institutions.

RECOMMENDED USE OF THE MEDICAL DEVICE as a part of complex therapy and rehabilitation after a stroke / acute disturbance of cerebral circulation;

after COVID-19;

subacute and chronic inflammation: prostatitis, benign prostatic hyperplasia, pneumonia, bronchitis, etc.;

poorly healing wounds and immune deficiency;

injuries and diseases of the joints and spine of various origins: arthritis, osteoarthritis, Para synovitis, epicondylitis, bursitis, back pain, sprains, bruises, myositis, tenosynovitis, etc.;

diseases of the cardiovascular system: primary and secondary hypertension, rheumatism and more;

diseases of the nervous system: splash, radicular syndrome, vibration disease, etc.;

inflammatory tissue diseases: mastitis, postoperative infiltration, strokes, etc.

Scope	Operating mode of the device	Recommended session duration	Recommended quantity of sessions
Modulation of the patient's immune system: restoration of the body's natural immunity.	Therapy	120 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Reduction of the rehabilitation period after surgery	Therapy / Rehabilitation	60 min in therapy mode. Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened. The next day 60 minutes in rehabilitation mode	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.

Reduction of pain (back, neck, rheumatic pain), reduction of spasmodic pain.	Therapy / Rehabilitation	60 min in therapy mode. Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened. The next day 60 minutes in rehabilitation mode	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Reduction of stress and nervous tension, prevention and prevention of colds and viral diseases	Therapy	20-30 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Reduction of stress and nervous tension, prevention and prevention of colds and viral diseases	Rehabilitation	40-60 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Rehabilitation after stroke	Rehabilitation	120 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Injuries and diseases of the joints and spine of various origins	Therapy	120 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.
Subacute and chronic inflammatory diseases, diseases of the cardiovascular system, diseases of the nervous system, diseases of the respiratory and auditory organs without purulent process, inflammatory skin diseases	Therapy	120 min Depending on the patient's condition, the time of therapy can be extended as prescribed by a doctor or shortened.	10 sessions in a row, one session daily, then a break of 7-10 days. A repeat course may be prescribed by a doctor.

1.3 CONTRAINDICATIONS



Proper examination and diagnosis must be performed, before starting treatment with the device. Please keep up to date with the latest developments and medical publications on devices with low-intensity electromagnetic radiation for detailed information on contraindications and side effects not known at the time of the device's manufacture. Contraindications listed in this section are given at the time of writing of the Instruction. No claims regarding the completeness of this list of contraindications are accepted. Before carrying out the procedures, a medical specialist should be convinced of the expediency of using this procedure, the responsibility for which he bears personally.

The use of the device is contraindicated if a patient has the following signs or pathologies:

- Bleeding, risk of bleeding and blood clotting disorders (haemophilia, haemorrhage, haemorrhoids and ulcers with the risk of bleeding, open wounds and injuries, etc.);
- Severe arterial obstruction (III and IV degree);
- Occlusive vascular diseases, such as obliterating arteriosclerosis and thromboangiitis obliterans (Buerger disease), in which organic occlusion and ischemia are detected;
- Swelling of tissues and the presence of foreign bodies in the affected area;
- Paroxysmal cardiac arrhythmia;
- Epilepsy;
- Gastric ulcer with a complicated course.
- The patient has implantable electronic devices (e.g, heart rate driver (pacemaker), etc.), as well as implantable devices containing metal parts.
- Individual intolerance to the procedure.

Pregnancy is an absolute contraindication to use. Individual intolerance of the procedure and/or discomfort during the procedure is an indication of its cancellation.

1.4 SIDE EFFECTS

Side effects have not been detected. Adverse effects of using the device are possible in case of neglect of the requirements for contraindications.

2 WARNINGS AND SAFETY PROVISIONS

2.1 SAFETY SIGNS ON THE DEVICE



The device is protected by reinforced insulation and has no galvanic connection to the ground

IP20

Degree of protection against external influences.

2.2 WARNINGS AND SAFETY PROVISIONS



The user must have the proper technical and medical qualifications and know the user's manual of this device in order to use this device. All maintenance procedures recommended by the manufacturer must be performed by personnel with appropriate approvals.



It is allowed to use the device in medical centers, in rehabilitation and sports medicine centers, SPA centers, massage rooms for adult patients (18 years and older).




The operator must inspect the housing of the electronic unit, as well as the power cord to ensure there is no external damage. Operation of the device with a damaged casing or a power cord is prohibited!



This device complies with the requirements of the electrical safety standard EN 60601-1: 2010.

It is necessary to connect the device to the mains supply in accordance with the national electrical safety regulations.



The device must be placed beyond the reach of a patient, especially children



Certified and safe materials are used for the device.



ATTENTION! Modification of the product is not allowed!



CAUTION! To avoid the risk of electric shock, the product must only be connected to a main supply that has a protective ground.



Disconnect the device from the power supply, before performing any cleaning or maintenance work. The means of simultaneous electrical separation of the supply circuits of the device from the circuits of the supply network is the mains switch of the device.



Connect the device only to a working socket with a rated voltage within the range 100 - 240V 50-60Hz. The location of the device should ensure that there is no tension on the power cord, unhindered connection and disconnection of the power cord from the mains are to be ensured to quick disconnection of the device from the mains in emergency situations.



Do not allow humidity to enter the electronics housing. Do not expose the device to dampness, vibration, or shock.



It is prohibited to use the device in a potentially explosive atmosphere, i.e. in the presence of a mixture of flammable anaesthetic gas with air, oxygen or nitrogen oxide. It is prohibited to use the device in rooms where flammable and potentially explosive substances are stored or used.



Potentially, there is a risk of passing microbes through the surface of the housing of the device. It is recommended to clean it regularly!



The patient should be properly located for treatment. It is necessary to monitor a patient's state during the procedure.



It is allowed to use the device only after studying this IFU!



It is prohibited to use the medical device in an oxygen-rich environment.

2.3 MEASURES TO PREVENT DAMAGE OF EQUIPMENT AND THE DEVICE

Connect the device to the mains through a circuit breaker with the characteristic "C" and a rated current of not more than 6A.

The placement of the device must ensure uninterrupted connection and disconnection of the power cord from the mains. Avoid the situation when the power cord is under the feet of a user or patient! Do not allow the mechanical load on the power cable and the device's enclosure (compression, stretching, stepping on, etc)!

It is prohibited to cover the device during operation.

It is forbidden to disconnect the device during operation from the mains network.

This device complies with the requirements of electrical safety and EMC (IEC 60601-1: 2005-12 3rd ed.). As a rule, the level of emitted electromagnetic interference is not sufficient to disrupt the operation of most devices. However, it should exclude the operation of the device in close proximity to sensitive equipment. It is recommended to place the device no closer than 3m to such equipment.

The device must be stored in a place protected from the direct sunshine.

It is necessary to exclude contact of the device with different solvents, gasoline, kerosene and other substances that can destroy or damage the device housing.

It is forbidden to install the device on slippery surfaces to prevent the device from falling.

3 BRIEF DESCRIPTION OF THE DEVICE

The therapeutic action on tissues and inner organs of a patient by low intensity electromagnetic field with frequency band of 100-1,500 MHz

Emitted electromagnetic waves result have the oscillatory effect in the human body, thus stimulating activity of the physical and chemical processes in the body. The penetrating ability of UHF waves in tissue is 8-11 cm on average. Skin and subcutaneous fat thickness have no significant influence on the reflection and absorption coefficients of the waves.

3.1 FUNCTIONAL DIAGRAM OF THE DEVICE

Studies of the effects of electromagnetic fields (EMF) on living organisms have been conducted since the middle of the last century. Thus, it was found that EMF can affect the biochemical reactions and behaviour of charged molecules near membranes, namely: to create electric fields in conductors, to exert force on moving charge carriers, to change the diffusion rate through membranes, to change valence angles, which affects binding proteins and synthesis of macromolecules, etc.

Studies in the field of molecular biology have established the presence of endogenous bioelectric signals, as well as determine their sources and effects on embryogenesis, regeneration and neoplasms. Ion fluxes and voltage gradients generated by ion channels and pumps are key regulators of cell proliferation, migration, and differentiation. Closed channels have movable folds in proteins, which in turn can be open, allowing ions to pass through the channel, or closed, preventing the passage of ions through the channel.

The uneven distribution of several key ions (Na^+ , Cl^- , K^+) between the intracellular and extracellular fluid and their movement across the plasma membrane determines the electrical properties of the membrane. All plasma membranes have a membrane potential, therefore, the membrane potential (V_{mem}) leads to the distribution of charges across the membrane. Each time the value of V_{mem} differs from 0 mV, in the positive or negative direction, the membrane is in a state of polarization. The magnitude of the polarization potential is directly proportional to the number of positive and negative charges separated by a membrane. In other words, changes in V_{mem} cause changes in the movement of ions across the membrane. Trigger events, such as the effect of an exogenous electromagnetic field (EMF), the frequencies of which resonate with endogenous EMFs, also cause changes in membrane permeability.

Changes in membrane potential regulate the proliferation of progenitor cells, stem cells and regenerative systems, as well as the efficiency of cytotoxic T lymphocytes.

Other studies have shown that stress gradients were not just membrane potentials but also specific signals for key metabolic processes in regenerative wound healing. These signals determine the path of cell migration, forming stress gradients between the intracellular and extracellular environment. Voltage gradients are localized DC electric fields that turn on and off at different stages of development. They spread to the extracellular space, as well as into the cytoplasm of one or more cells connected by slit compounds. These gradients can penetrate the cell membrane, the cytoplasm and even the membrane of the cell nucleus by signal transmission, with the EMF signal being received through receptors on the cell surface and then treated with G-proteins that bind receptors to effectors such as ion channels. It is known that these signaling processes have a correlation between the presence of EMF gradients and the cellular response.

In the process of technical implementation of the BIOL device, these features were taken into account and thus were determined before the application of the frequency range of meter and decimeter waves.


Functional diagram of the device BIOL (Pic. 3.1) includes radio frequency generator modules and module of emitters, operation mode processor, sensor of emission, control panel and power supply unit. Generator modules are designed to generate radio-frequency signals, amplify those signals and match them to the radiating aeriels, which are located in the radiator module. The signal is generated by frequency synthesizer microchip. The center radiated frequency value, deviation and the law of frequency variation are defined by the controller.

Control over the generator module controllers is carried through control bus of the operation mode processor. The switch-on sequence of the generators, session time, center frequency values and the law of their variation, as well as other parameters, are set according to the chosen operation mode. The control panel lets choose two possible operation modes: “Therapy” and “Rehabilitation”. The mode “Preparing” is run automatically after each switching on the device and provides self-testing and diagnostics of the hardware.

Using the control panel, you can switch on/off the apparatus, choose operation mode (Therapy/Rehabilitation), and set the session time. The control panel indicator is used to visualize parameters of the assigned mode of operation, and, after the beginning of the session, time to finish the session and the total time of the device’s working.

A built-in sensor of emission provides sequential control over the emitted radio signal in the “Preparing” mode.

A power supply unit is designed to generate stabilized DC voltage of 5 V. Voltage converters (V/C) of the generator modules provide the formation of boosting DC voltage, which is necessary for the operation of the modules. Radio frequency amplifiers (RFA) are designed to amplify a signal, generated by the signal synthesizers, to the rated power value.



Picture 3.1 Functional Diagram of the device BIOL

4 BEFORE USING THE DEVICE BIOL

4.1 THE LIST OF PREPARATORY ACTIONS

Before using the device BIOL, it is necessary to do the following:

- ✓ Make sure that the power supply network is used with a protective ground; that the voltage in the network is in the range 100 - 240V, 50 - 60 Hz. The device BIOL is intended for connection to type F (Schuko) sockets (European socket with CEE 7/4 grounding, DIN 49440 standard). A suitable adapter is required to connect the device to other types of sockets; in any case, the presence of a protective ground in the power socket is mandatory!
- ✓ Ensure the presence and storage of 70 - 96% water-alcohol solution to clean the device.
- ✓ Ensure the presence and storage of medical alcohol wipes or cotton pads for cleaning the device.

- ✓ Ensure the presence and storage of wet wipes for a screen that does not contain alcohol to clean the display of the device BIOL.
- ✓ Organize the workplace of the operator so that the device is placed on a solid, smooth, dry and not slippery surface to comply with the measures listed in paragraph 2.3.
- ✓ Remove the device from its packaging. Check if there is no damage to the housing of the device and the power cord. Check that the power switch is in the position "Off" (the button on the front panel is in the "not pressed" position.)
- ✓ Connect the power cord to the power connector on a back panel of the electronic unit, plug the power cord into a power socket.

4.2 THE OPERATOR'S QUALIFICATION

The device BIOL is intended for use by the operators having special knowledge in the field of application of this device and trained about the proper application of the device, as well as the operators who have practical skills in working with similar medical equipment.

The operator should have main physical and cognitive abilities, such as sight, hearing and literacy. A tremor in hands of the operator is an obstacle for the device use, as the parameters of a session could not be set.

Besides, the operator must take into account the manufacturer's recommendations (Chapter 1, item 1.2 "Indications for use" and item 1.3 "Contraindications") to be aware of the latest developments and medical publications for detailed information on contraindications and side effects, not known at the time of manufacture.

The operator has to take the appropriate training regarding the correct operation of the device before working with it:

- ✓ Intended use of the device with practical exercises;
- ✓ The mechanism of action and function of the device;
- ✓ Setting up of the working modes;
- ✓ Recommendations for use of the device;
- ✓ Contraindications and side effects;
- ✓ An explanation of alerts in all modes of operation;
- ✓ Method of functional verification of the device.

Further recommendations for the scope of training may vary depending on the country. Please contact your local representative of SC KOLIBRI LLC for detailed information on training.


5 OPERATION OF THE DEVICE BIOL



In a process of the device application, the operator must be in a satisfactory physical and emotional state (after a sufficient rest), should not take in psychotropic substances, analgesics, opiates, sleeping pills, drugs or alcohol for at least 48 hours before performing procedures using the device.

5.1 A PATIENT'S LOCATION

The patient should be at a distance of 1.5-2 meters from the device from any side in a convenient position (sitting, lying, etc.).



5.2 AN OPERATOR'S LOCATION

The operator is not recommended to be within a radius of 3m around the device during the treatment session.

5.3 RECOMMENDED OPERATING MODES

The device BIOL is pre-programmed for two main working modes, namely:

«THERAPY»	«REHABILITATION»
The mode, which is used for treatment	The mode, which is used for prevention of diseases

5.4 SWITCHING THE DEVICE ON

- 5.4.1 Connect the power cord to the device and plug it into the power supply.
- 5.4.2 Press the power button on the front panel of the device:



The red indicator on the button will light up


- 5.4.3 The information on the process of preparing the device for work will be displayed ("Preparing"), during which the automatic self-testing and diagnostics of the device are carried out:



5.4.4 Upon completion of preparing, the device is ready to be used "Ready" and it goes into the standby mode of setting the session parameters. The information about the total operating time of the device is displayed "Total: XXh XXm" (total operating time in hours (h) and minutes (m) since it was manufactured):



When preparing the device for work, during the self-testing and diagnostics ("Preparing" mode), 8 rectangles must be displayed successively.



In the absence of at least one of the rectangles, the device is to be repaired. Further use of the device is prohibited. Contact the manufacturer!

5.5 DESCRIPTION OF CONTROL MEANS

The only one control element of the device is the knob (encoder), which rotates both clockwise and counter clockwise (without restriction), and also pressed (as a button). It is located on the front of the device to the right of the display.



5.6 SETTING UP THE PARAMETERS OF SESSION



To prevent damage to the knob-encoder, do not apply excessive force when pressed, avoid impacts and strokes. Keep the device in the supplied bag. Put a damper compactor around the knob-encoder when transporting (included in the delivery set).

5.6.1 Setting the operating mode

Press the knob-encoder perpendicular to the front panel, the same as you press the normal button. The Menu will appear on the display: «Rehabilitation» or «Therapy» mode will be displayed



To switch the mode from “Therapy” to “Rehabilitation” or vice versa, roll the knob-encoder left or right. Once the menu you need is displayed, press the encoder-knob, as a normal button. The operating mode of the device is set.

5.6.2 Setting the time of session

After setting up the mode of operation, the menu for time setting will be displayed (duration of a session).



First you need to set the duration of the session in **hours** - «Enter hours...». To start, press the encoder-knob, as a normal button. Then rotate the knob-encoder left or right to select the required number of hours for a session. Set hours by pressing the knob-encoder, as a button.

Notice: If the duration of the session does not exceed one hour, then set the number 0 (zero) to the "Enter hours ..." field, according to the instructions above.

The next step is to set the number of **minutes** of the required operating time – it is done the same as for the hours, by rotating and pressing the knob-encoder (described above).

After setting the time of operation (both hours and minutes), the session begins and the information about set mode ("Therapy" or "Rehabilitation") and time of session "Time: " is displayed.



The operating time of the device on the display will gradually decrease (the timer works in the opposite direction): from the set number to zero.

At the end of the set working time (session), a short beep will sound and the unit will go into the "Ready" mode automatically. The device is ready again to set the parameters of a new session.

5.7 SHUTDOWN AND STORAGE

To shut down the device BIOL properly, press the power button, which is located on the front panel of the device to the right. The LED on the button goes out. Then you can unplug the power cord and place the device in a portable protective case and put a damper compactor around the knob-encoder (included in the delivery set).



5.8 MESSAGES ON THE DISPLAY



Automatic mode of preparation of the device to work, which includes self-testing and diagnostics of hardware are carried out – “Preparing”



Standby mode - the device is ready for setting the mode of operation and time of session – “Ready”




The mode for diseases prevention - «Rehabilitation»



The mode for diseases treatment - «Therapy»



Menu to setup the time of a session



Indication of the device's operation in the set mode «Rehabilitation»



Power indication of the device (the red indicator on the button is light when the device is turned on)..



Power indication of the device (the red indicator on the button isn't light when the device is turned off).

Sound signal

The short beep is sound when any of the operating mode is completed and the device is switched to the “Ready” mode.

6 MAINTENANCE OF THE DEVICE

6.1 CLEANING

Regular cleaning of the device ensures its reliable and trouble-free operation.

Before cleaning and/or repair, disconnect the device from the mains.

In general, external cleaning of the device's housing is carried out depending on the frequency of the device use.

All the details that are in contact with the operator must be cleaned with medical alcohol wipes or cotton pads soaked in 70-96% water-alcohol solution.

It is very important to avoid an enter of liquids inside the device.

It is necessary to keep clean the ventilation slots of the device.

Only special non-woven wipes for LCD monitors are allowed to clean the TFT display. Those wipes should not contain alcohol.

6.2 FUSE REPLACEMENT

If the device does not work, after connecting to a power supply and turning on, so it is necessary to check and replace a fuse/fuses: maybe one of them is failed (swelled). The fuse holders are located on the back panel of the device. Please follow these steps to replace the fuse:

- ✓ Unplug the power cord from the power socket.
- ✓ Unlock one of the fuse holders on the panel.
- ✓ Remove the damaged fuse from the holder.
- ✓ Install a new fuse in the holder (type T4AL250V or T4AH250V).
- ✓ Lock the fuse holder back into the panel



6.3 MAINTENANCE AND SAFETY CHECK

Preventive maintenance is not necessary. However, regular maintenance can help identify possible defects at an early stage and thus increase safety and extend the life of the device.

It is recommended to perform functional checks and safety checks of the device at least once a year. If the national safety regulations for medical devices require a more frequent periodicity of

tests and inspections, it is necessary to follow national regulatory documents. Functional and safety checks are carried out at the producing factory or authorized service centres.

6.4 DISPOSAL AND ENVIRONMENTAL PROTECTION



In case of failure and impossibility of further use of the device BIOL, it is disposed of as used electronic equipment. Please dispose of the apparatus in accordance with the current regulations in your country.

PACKAGING DISPOSAL: Packaging components (cardboard, expanded polystyrene, etc.) are classified as solid waste and therefore they can be easily recycled by using recycling processes. Before sending the components to special recycling centers, it is necessary to check local regulations in that regard in order to comply with them fully.

PRODUCT DISPOSAL: The device BIOL consists of various materials. Nevertheless, all of them (metal, plastic, electrical conductors, printed circuit boards, chips, etc.) do not contain hazardous substances and they can be sent to special recycling centers in the same way as electronic equipment. Before sending the components to special recycling centers, it is necessary to check local regulations in that regard in order to comply with them fully.

6.5 REPAIR

Only personnel who has an appropriate authorisation from the company Scientific company KOLIBRI LLC can carry out repairs of faulty devices BIOL. Only the spare parts indicated by Scientific company KOLIBRI LLC must be used for that. The personnel with appropriate authorisation may include Scientific company KOLIBRI LLC staff and technician specialists of sales representatives who have permission from Scientific company KOLIBRI LLC.

6.6 LIFETIME

Considering the characteristics of similar equipment on the market, as well as the actual period of the device BIOL being marketed (since 2013), the following lifetime period is established for the device, with the obligatory compliance with conditions of packaging, storage, transportation and use:

- 5 years or 2 500 working hours - for the electronic unit of the BIOL;

The probability of failure of the components and accessories of the device increases after exceeding of the lifetime.

7 COMPLETION AND RECOMMENDED SERVICE MATERIALS

7.1 COMPLETE SET OF THE DEVICE

The device BIOL complete set includes the following:

- ✓ Electronic unit BIOL;
- ✓ Power cable;
- ✓ Bag (case) for transportation and storage of the device BIOL;
- ✓ Instruction for Use with warranty card

7.2 RECOMMENDED MATERIALS FOR SERVICE AND MAINTENANCE

It is recommended to use the following items for cleaning:

- ✓ 70 - 96% water-alcohol solution and cotton pads;
- OR
- ✓ Medical alcohol wipes

It is recommended to use wet wipes for monitors that do not contain alcohol to clean the display of the device BIOL.

8 DEVICE APPEARANCE



9 TECHNICAL DESCRIPTIONS

9.1 CLASSIFICATION

By the way of protection against electric shock, the device belongs to the class I with a working part of type BF.

The degree of protection against external influences is IP20.

The operating mode of the electronic unit of the device is long.

In terms of the degree of risk of medical use, the device belongs to class IIa according to the Technical Regulation No. 753 and the Directive 2007/47/ EC which amends the Directive 93/42 / EEC.

9.2 TECHNICAL DATA

The electronic unit	
Input voltage of the network	100-240 VAC 0.35-0.2 A
Network frequency	50/60 Hz
The mains fuse type	T4AH250V or T4AL250V
Electricity consumption	max. 15V.A
Total output power	0.1 W
The ambient temperature during the operation process	5 – 40 °C
The ambient temperature during a storage and transportation	-25 ° C without relative humidity control + 70 ° C with relative humidity control This class 7K3 as described in IEC/TR 60721-4-7:2001
Ambient air pressure	700-1060 hectopascal
Air humidity	15%-93%, without condensation

Total weight	2,0 kg.
Net weight	1,5 kg.
Dimensions of the device BIOL (D / W / H)	260x180x65 mm.
Protection against water penetration	IP20
Software version	V.1.0

The mains switch of the device is the means of simultaneous electrical separation of the power circuits of the device BIOL from the supply network circuits.

The temperature of all surfaces of the medical device does not exceed $40 \pm 0.5^{\circ}\text{C}$

9.3 TRANSPORTATION AND STORAGE

Transportation and storage of the device BIOL are only permitted in the manufacturer's packaging.

You should avoid shaking and impacting of the package during transportation and storage.

Storage and transportation conditions of the device are the following:

- ✓ from -25°C (without relative humidity control) up to $+70^{\circ}\text{C}$ (with relative humidity control);
- ✓ relative humidity 15% -93% without condensation;
- ✓ And also, the absence of aggressive impurities that cause corrosion in the air.
- ✓ Storage 7 years.

10 WARRANTY TERMS

The warranty for the device BIOL is 36 months from the date of sale.

The warranty does not cover cables and fuses.

The manufacturer or the authorized representative provides the free repair of malfunctions or replacement of the device during the warranty period if there is a detection of manufacturing defects or defects in materials. The warranty for such a device does not apply if faults are caused by the user, due to a violation of the operating rules of the device described in this IFU, or the device was misused.

The warranty also does not cover damages caused by the user violating the storage and transportation rules set out in this IFU, as well as the force majeure.

Warranty claims are accepted only on the condition that the device is returned in its complete configuration, in its pure form, without external mechanical damages and traces of disassembly/opening.

Transportation costs and the risk of accidental loss during the delivery of the returned product shall be borne by the customer.




ATTENTION! It is not allowed to make any changes in the design of the device. Any unauthorized opening, repair or modification of the system by unauthorized personnel releases the manufacturer from obligations and responsibility for the safe operation of the device. In this case, the warranty is automatically declared invalid even before the expiration of the warranty period. The warranty is cancelled if the customer has made a modification or made any uncoordinated changes to the software of the device without the written consent of the company Scientific company KOLIBRI LLC.

For all questions regarding the operation of the device, please contact our customer service:
Scientific company KOLIBRI LLC

Provulok Kinnyi 8a, 61001 Kharkiv, Ukraine
email: info@kolibri.one
tel. +380913011110

11 POSSIBLE ERRORS AND MALFUNCTIONS. TROUBLESHOOTING

	Error / problem	Possible reasons	Indicators / Disposal
E1	The device does not turn on	No power supply	-Check the presence of an electrical current in the mains. In case of absence, try again later. - Check the connection of the power cord to the device and to an electrical outlet. - Check the network cable for defects if there are any and replace it.
		The refractory fuses are defective or missing	Check the installed refractory fuses; replace the damaged ones or install the missing ones.
E2	Less than 8 indicator rectangles appeared on the display after run the device, in the process of Preparing 	One (or more) internal generator is defective.	The device must be repaired at the factory. Send the device to the nearest authorized service center or representative office of SC KOLIBRI LLC or to the manufacturer.

12 CHECKLIST

No	Date of verification / control Warranty / post-warranty service	Remark	Note	Performer	Sign
1.					
2.					
3.					

13 PRODUCT ACCEPTANCE PROTOCOL

The device BIOL, No **B** _____ complies with the technical requirements and is completely serviceable.

The warranty period is 36 months from the date of delivery of the equipment.



Date of manufacture “ _____ ” _____ 20____.


(signature)

(First name and Last name)



SC KOLIBRI LLC
Provulok Kinnyi 8a, Kharkiv,
61001, Ukraine

Tel.: +380913011110
email: info@kolibri.one
URL: <https://kolibri.cloud>



Authorized representative in European Union: “ONKOCET Ltd.”

4 Kutuzovova str., 90201 Pezinok, Slovakia,
tel.: +421 (2) 44 64 09 77
email: onkocet@onkocet.eu
URL: www.onkocet.eu

14 PACKING LIST

The device BIOL, No **B** _____ is packed in the company Scientific company KOLIBRI LLC in accordance with the technical requirements.

Date of packing « _____ » _____ 20____.

A packer:


Signature

First name and Last name



SC KOLIBRI LLC
Provulok Kinnyi 8a, Kharkiv,
61001, Ukraine

Tel.: +380913011110
email: info@kolibri.one
URL: <https://kolibri.cloud>




Authorized representative in European Union: “ONKOCET Ltd.”

4 Kutuzovova str., 90201 Pezinok, Slovakia,
tel.: +421 (2) 44 64 09 77
email: onkocet@onkocet.eu
URL: www.onkocet.eu

15 WARRANTY CERTIFICATE

For repair (replacement) during the warranty period
Medical Equipment - **BIOL**

Serial number of the device B_____		Date of manufacture «_____» _____ 20__
Date of purchase «_____» _____ 20__		Signature and seal of the seller
Date of commissioning «_____» _____ 20__		Signature
The operating time of the device _____ [hour:min]		Signature

The warranty for the device BIOL is 36 months from the date of sale.

The warranty does not cover power cables and fuses.

The manufacturer or his authorized representative performs free repair of malfunctions or replacement of the device during the warranty period, if there is detection of manufacturing defects or defects in materials. The warranty for such a device does not apply if there are faults caused by the user, due to a violation of the operating rules of the device described in this manual, or the usage of the device for other purposes.


The warranty also does not cover damages caused by the user violating the storage and transportation rules set out in this manual, as well as the force majeure.

Warranty claims are accepted only on condition that the device is returned in its complete configuration, in its pure form, without external mechanical damages and traces of disassembly / opening.

Transportation costs and the risk of accidental loss during the delivery of the returned product shall be borne by the customer.




SC KOLIBRI LLC
Provulok Kinnyi 8a, Kharkiv,
61001, Ukraine
Tel.: +380913011110
email: info@kolibri.one
URL: <https://kolibri.cloud>






Authorized representative in European Union: "ONKOCET Ltd."
4 Kutuzovova str., 90201 Pezinok, Slovakia,
tel.: +421 (2) 44 64 09 77
email: onkocet@onkocet.eu
URL: www.onkocet.eu

All valid permits and certificates are available on the manufacturer's website <https://kolibri.cloud/>

16 LABELLING



17 THE COMPLETE SET OF BIOL

Name of a part	Quantity
	The unit of MD BIOL 1
	Power cable (Type C13) 1
	Package 1
IFU (offered in electronic form/available on web)	1

18 DOCUMENT HISTORY AND VERSION CONTROL

Version	Version Date	Summary of changes	Author	Related documents
1.0	2021-01-12	Created	Team, listed on the title	



Read the IFU carefully before using the device!